



Accessories (each sold separately):

- Clip-On Caddy Dispensers: MP-015 (below left) and MP-085
- BOWMAN® Door and Wall Hangers: HH004-0412, MB-450, MB-500, MB-510, and MB-600 (below right)



Product Description:

- · Protection System
- Holds personal protection equipment in two bulk compartments, one box of earloop masks, shields, or boxed gowns, disposable or reusable gowns, aprons, or coveralls/jumpsuits, and three boxes of gloves
- · Side load bulk and gown areas
- Keyholes for wall mounting, or door hanger cut out for use with wall or door hanger
- Quartz Beige ABS Plastic

Primary Area for Product Usage:

• Anywhere personal protection equipment is needed

Dispenses:

- Top left and right: Item(s) no larger than
 7.90"W x 10.00"H x 4.55"D (20.1 cm x 25.4 cm x 11.4 cm)
- Top center: Box(es) no larger than 8.05"W x 10.00"H x 4.55"D (20.4 cm x 25.4 cm x 11.6 cm)
- Bottom left: Item(s) no larger than 14.15"W x 16.33"H x 4.55"D (35.9 cm x 41.5 cm x 11.6 cm)
- Right: Box(es) no larger than 10.00"W x 5.38"H x 4.55"D (25.4 cm x 13.7 cm x 11.6 cm)


Package Specifications:

- Package Quantity: 1 per case
- 26"L x 6"W x 32"H (66.0 cm x 15.2 cm x 81.3 cm)
- 11.0 lbs (5.0 kg) approximated

Product Specifications (overall external dimensions):

- 24.47"W x 26.63"H x 4.75"D (62.2 cm x 67.6 cm x 12.1 cm)
- 8.2 lbs (3.7 kg) approximated

Mounting Specifications:





PS034-0212 Protection System

Installation Instructions/Characteristics & Care

Installation

Always mount dispensers securely prior to loading with product.





Parts included:

Dispenser

Tools Required:

· No tools required to assemble dispenser

Care and Construction ABS Plastic

Material Characteristics

Designed for Long-Term, Heavy-Duty Applications

The combination of copolymers (Acrylonitrile, Butadiene, and Styrene) gives ABS an excellent combination of strength, rigidity, and toughness. The fauxwood finish is ARVINYL laminate created from a semi-rigid film.

ABS is chemical and stress-crack resistant. It is tough, dimensionally stable, and will hold up to long-term use. To ensure product longevity, fauxwood products should not be mounted in very hot or cold areas.

Care and Cleaning

Use mild soap and water to clean. Petroleum-based oils, solvents, and paints should not be used.

Standard hospital cleaning supplies are approved for use on our ABS plastic products.